


MOTECK®




醫療病床 / 居家照護床用

Product Catalog for Medical Bed / Home Care Bed


產品組目錄



主控制盒
Control Box



搭配推桿
Actuator



控制器
Control





三功能床 *3 Functions Bed*



● 居家照護

Home Care

主控制盒 Control Box



CBT2-230VAC-30VDC-3

- 輸入電壓 Input voltage
100 / 120 / 230 V AC (transformer) /
100 ~ 240V AC (SMPS)
- 輸出電壓 Output voltage
24 V DC (standard)
- 驅動推桿數 Number of actuator
1 / 2 / 3
- 殼體顏色 Color
黑 Black 
淺灰 Light gray RAL 7001 

搭配推桿 Actuator



FD60



背部升降 Back
FD60-24V-A4-370.570



腿部升降 Leg
FD60-24V-A4-273.373



床體升降 Hi-Low
FD60-24V-A4-323.473

控制器 Control



HC8-AD-B2

配件 Optional Accessory



可插拔電源線
Exchangeable cable



三功能床

3 Functions Bed

● 醫療病床

Medical Bed


主控制盒 Control Box



MD6C-230VAC-24VDC-1-3

- 輸入電壓 Input voltage
100 / 120 / 230 V AC
- 輸出電壓 Output voltage
24 V DC
- 驅動控制器數 Number of control
1 / 2
- 驅動推桿數 Number of actuator
1 / 2 / 3 / 4

搭配推桿 Actuator



MD66



背部升降 Back
MD66-24-A4-370.570




腿部升降 Leg
MD66-24-A4-273.373



床體升降 Hi-Low
MD66-24-A4-323.473

控制器 Control



H3T-08-3

配件

Optional Accessory



可插拔電源線
Exchangeable cable



附控制面板三功能床

3 Functions Bed with control panel

● 醫療病床

Medical Bed

主控制盒 Control Box



MD6C-M-SMPS-24VDC-1-4

- 輸入電壓 Input voltage
100~240V AC
- 輸出電壓 Output voltage
24 V DC
- 驅動控制器數 Number of control
1
- 驅動推桿數 Number of actuator
2 / 3 / 4

搭配推桿 Actuator



MD60



背部升降 Back
MD60-24V-A8-370.570



MD66



腿部升降 Leg
MD66-24-A4-273.373



床體升降 Hi-Low
MD66-24-A4-323.473

控制器 Control



CP1-20-4



床尾板控制器
Foot Control



床側板控制器
Side Rails Control

配件

Optional Accessory



JB1集線器 (床底燈)
Junction Box with or
without Under Bed Light



BTP2充電電池盒
Battery Box



五功能床 *5 Functions Bed*

● 居家照護

Home Care

主控制盒 Control Box



CB4M-230VAC-24VDC-1-4

- 輸入電壓 Input voltage
110 / 230 V AC
- 輸出電壓 Output voltage
24 V DC
- 驅動控制器數 Number of control
1 / 2
- 驅動推桿數 Number of actuator
2 / 3 / 4

搭配推桿 Actuator



FD60



背腿同動 Joint Movement



背部升降 Back

FD60-24V-A4-370.570




腿部升降 Leg

FD60-24V-A4-273.373



床體升降 Hi-Low




前傾 / 後傾 Forward-Backward

FD60-24V-A4-323.473

FD60-24V-A4-323.473

控制器 Control



HC10-AE-B2

配件

Optional Accessory



可插拔電源線
Exchangeable cable




五功能床

5 Functions Bed

● 醫療病床

Medical Bed


主控制盒 Control Box








CB4P-230VAC-24VDC-1-4

- 輸入電壓 Input voltage
110 / 230 V AC
- 輸出電壓 Output voltage
24 V DC
- 驅動控制器數 Number of control
1 / 2
- 驅動推桿數 Number of actuator
3 / 4


搭配推桿 Actuator



MD66

- 背腿同動 Joint Movement

MD66-24-A4-370.570
- 背部升降 Back

MD66-24-A4-273.373
- 腿部升降 Leg

MD66-24-A4-273.373
- 床體升降 Hi-Low

MD66-24-A4-323.473
- 前傾 / 後傾 Forward-Backward

MD66-24-A4-323.473

控制器 Control



H3T-10-4

配件 Optional Accessory



可插拔電源線
Exchangeable cable




附控制面板五功能床

5 Functions Bed with control panel

● 醫療病床

Medical Bed







帶快速釋放五功能床 *5 Functions Bed with Quick Release*

● 醫療病床

Medical Bed





翻身床 *Turning Frame Bed*

● 居家照護

Home Care

● 醫療病床

Medical Bed


主控制盒
Control Box



CB4M-230VAC-24VDC-1-4

- 輸入電壓 Input voltage
110 / 230 V AC
- 輸出電壓 Output voltage
24 V DC
- 驅動控制器數 Number of control
1 / 2
- 驅動推桿數 Number of actuator
2 / 3 / 4

搭配推桿
Actuator



FD60



背部升降 Back

FD60-24-A4-373.538



腿部升降 Leg

FD60-24-A4-323.448



床體左右翻 Left-Right


FD60-24-A4-365.535



便盆 Bedpan

FD60-24-AG-323.488

控制器
Control



HC-10-4



多功能秤重床

Multiple Functions Bed with Load Cell Systems

● 醫療病床

Medical Bed


主控制盒 Control Box



CB5P-M-230-24

- 輸入電壓 Input voltage
100 / 120 / 230 V AC
- 輸出電壓 Output voltage
24 V DC
- 驅動控制器數 Number of control
2
- 驅動推桿數 Number of actuator
5

搭配推桿 Actuator



MD60+QR1



背部升降 Back
(帶快速釋放 with Quick Release)
MD60-370.570


控制器 Control



CP2



H3T-M




床側板控制器
Side Rails Control

轉接盒 Junction Box



JB3
Junction Box



MD62



腿部升降 Leg
MD62-273.373




床體升降 Hi-Low
MD62-323.473
MD62-323.473

配件 Optional Accessory



FL2床底燈
Under Bed Light



TLC1秤重系統
Load Cell Systems

LOCATION

總公司 Head Office

台灣新北市汐止區新台五路一段79號1樓之1
1F-1, No. 79, Sec. 1, Xintai 5th Rd., Xizhi Dist.,
New Taipei City 22101, Taiwan
Tel : +886-2-2698-1220
Fax : +886-2-2698-1324



吳江工廠 Suzhou

中國江蘇省吳江市汾湖經濟開發區湖北路500號
No.500, Hubei Road, Fenhu Economic Development Zone,
(215213) Wujiang City, Jiangsu Province, China
Tel : +86-512-6329-6328
Fax : +86-512-6329-6325

廣東辦事處 Gangzhou

中國廣東省東莞市寮步鎮石大路銀基大廈503室
Room 503, Yinji Building, Shida Road, Liaobu Town,
Dongguan, Guangdong, China
Tel : +86-769-8375-9338
Fax : +86-769-8375-9328

