

Bed Side Rail BP21, BP22

BP21, BP22 Bed Side Rail can be used as an extension of the control system. It helps nurses and patients to easily operate first aid, care, bed lifting and other functions, and is suitable for medical institutions and home care.



Features and Option


Main applications: Medical care

- Available number of buttons
 - BP21 side rail (patient side): 4
 - BP21 side rail (nurse side): 12
 - BP22 foot board: 15
- Weight
 - BP21 side rail (back): 2.6±0.2 KG
 - BP21 side rail (leg): 2.2±0.2 KG
 - BP22 head board: 3.8±0.2 KG
 - BP22 foot board: 4.0±0.2 KG
- BP21 bed lifting angle indication
- Power on indicator
- Bed low position warning indicator
- With 4 red angle beads
- Material: PP
- Color: Pantone 9100 C
- Customized panel


Dimensions

Unit: mm


- BP21 side rail (back)




- BP21 side rail (leg)



- BP21 bracket





- BP22 head board & foot board





Button Definition

- BP21 side rail (patient side)




No.	Definition
①	Back Up
②	Back Down
③	Leg Up
④	Leg Down

- BP21 side rail (nurse side)




No.	Definition	No.	Definition
①	Lock Button	⑦	Height adjustment: Up
②	CPR	⑧	Height adjustment: Down
③	Lateral Tilt Left	⑨	Leg Up
④	Lateral Tilt Right	⑩	Leg Down
⑤	Reverse Trendelenburg	⑪	Back & Leg: Up
⑥	Trendelenburg	⑫	Back & Leg: Down

- BP22 foot board



Lock button indicator (orange)



No.	Definition	No.	Definition
①	Under Bed Light	⑨	Back & Leg: Up
②	CPR	⑩	Back & Leg: Down
③	Lateral Tilt Left	⑪	Leg Up
④	Lateral Tilt Right	⑫	Leg Down
⑤	Reverse Trendelenburg	⑬	Back Up
⑥	Trendelenburg	⑭	Back Down
⑦	Height adjustment: Up	⑮	Lock Button
⑧	Height adjustment: Down		

Standard Panel Artwork


- Code A (with Patient Turnover function)

Parts	MOTECK P/N	Artwork
BP21 side rail (nurse side)	52H0740	
	52H0742	
BP21 side rail (patient side)	52H0741	
	52H0743	
BP21 side rail (leg)	52H0755	
	52H0756	
BP22 head board	52H0757	
BP22 foot board	52H0744	

- Code B (without Patient Turnover function)


Parts	MOTECK P/N	Artwork
BP21 side rail (nurse side)	52H0777	
	52H0778	
BP21 side rail (patient side)	52H0741	
	52H0743	
BP21 side rail (leg)	52H0755	
	52H0756	
BP22 head board	52H0757	
BP22 foot board	52H0781	

System Diagram




Cable Length Diagram

- BP21




- BP22




Installation


- Bed frame assembly diagram




- BP21 Up



- BP21 Down



- BP22 installation diagram



Ordering Key

[illegible]

Terms of Use

The user is responsible for the suitability of MOTECK products, and the products listed on the MOTECK website are subject to change without notice. MOTECK reserves the right to terminate sales or delete any products displayed on the website or listed in its catalog.